


TBS INOX[®] USA

DELIVERING PREMIUM PRODUCTS
"FOR BETTER DENTISTRY"
SINCE 1992



ABOUT US



QUALITY INSTRUMENTS

MISSION STATEMENT

Deliver premium, innovative products "For Better Dentistry" using our global resources to set standards in design, quality, and efficiency.

COMPANY OVERVIEW

TBS Instruments originated 30 years ago, with merely a small team selling quality dental products at an affordable price to dentists across North America. Over the years, with its vast manufacturing and distribution experience, the company has transformed into a powerhouse of innovation and excellence.

With thorough research and development and an intent of modernization under way, TBS has designed a number of high quality dental instruments, with its prominent invention being FRINGS – forceps with spring. Being the first and only company in the world to have formulated forceps with auto retractable technology, TBS has set a new standard of revolutionary instruments with optimal performance in the dental industry.

Today, TBS has reached dentists internationally, garnering KOL's and a clientele throughout the world. With its team now expanded to 500+ individuals, TBS hopes to continue on the path of advancement, while creating superior and groundbreaking dental instruments.

KEY OPINION LEADERS

Dr. Neal Patel
Colombus, Ohio



Dr. Isaac D. Tawil
Brooklyn, N.Y.



Dr. Scott Ganz
Fort Lee, New Jersey



Dr. Bart Silverman
New City, N.Y.



Dr. Miguel Stanley
Lisbon, Portugal



Dr. Achraf Souayah
Tunisia



Dr. Irfan Abas
Amsterdam, Holland



Dr. Mitsias Miltiadis
Athens, Greece

FRINGS®

forceps with spring

FRINGS® feature a never before seen, spring-driven "auto-retractable" design, so you no longer have to use your fingers to open the forceps. EVER. They are available in a wide range of the most desirable beak designs, for optimal grip – reducing the risk of tooth fracturing. FRINGS® are crafted from 100% German stainless steel and feature a trademarked two-tone blue finish with an ergonomical handle that will change your dental practice forever! And, FRINGS® are protected by a lifetime warranty.



Visit our website
tbsinox.com
or call us
+1-844-TBS-INOX (827-4669)
to learn more!

FOLLOW US



/TBSINSTRUMENTS

TBS INOX®
USA

ERGONOMIC
LIGHTWEIGHT
HANDLE
provides
better grip
and control.

LIFE TIME
WARRANTY

100% CRAFTED WITH
**GERMAN
STEEL**

EXCLUSIVE
TRADEMARKED
2 TONE BLUE
FINISH



AUTO-RETRACTABLE





STANDARD

33001F

ANTERIOR,
TAPERED NOTCHED
SERRATED BEAKS



150

33002F

UPPER UNIVERSAL,
TAPERED
SERRATED BEAKS



79

33003F

LOWER MOLARS,
BEVELED, TAPERED
SERRATED BEAKS



APICAL

33004F

UPPER PREMOLARS
& CANINES, BEVELED,
TAPERED SERRATED BEAKS



67

33007F

UPPER MOLARS,
TAPERED SERRATED
BEAKS



**23
COWHORN**

33008F


LOWER 3RD MOLARS,
SPECIAL TAPERED
HORNS



**45
ROOTS**

33009F

LOWER ROOTS,
TAPERED SERRATED
CUPPED BEAKS



51 ROOTS

33010F

UPPER ROOTS,
TAPERED
SERRATED CUPPED BEAKS

FRINGS®
forceps with spring

**SPECIAL
 EDITION**



33005F
TAWIL
 MODIFIED 151
 LOWER UNIVERSAL,
 TAPERED, BEVELED,
 SERRATED BEAKS



33006F
PATEL
 MODIFIED COWHORN
 LOWER THIRD
 MOLARS, TAPERED,
 SERRATED COWHORNS



33011F
GANZ
 MICRO
 ROOTS,
 CROSS SERRATED BEAKS



Special fine serrated, titanium curettes designed for curettage, granular tissue removal and tooth socket debridement are optimal for any surgical procedure. TBS serrated curettes are loved by clinicians all around the world. Our serrated curettes are available in four different sizes.

PERIODONTAL SERRATED SURGICAL CURETTES



LARGE 12 4mm WIDE SERRATED SPOON **38001**



SMALL 10 3mm WIDE SERRATED SPOON **38002**



MEDIUM 11 3.5mm WIDE SERRATED SPOON **38003**



MINI 9 2.5mm WIDE SERRATED SPOON **38004**



FINE SERRATION on the blade is optimal to cut the periodontal ligament.



SHARK HANDLE™ exclusive design to provide optimal grip

100% CRAFTED WITH
**GERMAN
STEEL**

PERIODONTAL ELVATOMES

Elvatomes™ have the sharpness and precision of periostomes and sturdiness of elevators. Spade fine and very sharp point are optimal to enter skintight spaces and cut the periodontal ligament for an atraumatic extraction.



38021 STRAIGHT



38022 DISTALLY ANGLED



38023 MESIALLY ANGLED



38024 LEFT SPADE



38025 RIGHT SPADE



38100 ELVATOME™
KIT (5 PACK)

Visit our website
tbsinox.com
or call us
+1-844-TBS-INOX (827-4669)
to learn more!

FOLLOW US



/TBSINSTRUMENTS

TBS INOX
USA

100% CRAFTED WITH
GERMAN
STEEL

PERIODONTAL PERIOTOMES



FLEXIBLE, MICRO SERRATED,
TITANIUM NITRIDE

STRAIGHT 38017



FLEXIBLE, MICRO SERRATED,
TITANIUM NITRIDE

ANGLED 38018




FLEXIBLE, MICRO SERRATED,
TITANIUM NITRIDE

REVERSED 38019



NON-FLEXIBLE, MICRO
SERRATED, TITANIUM NITRIDE



CONCAVED 38020



**PERIOTOME
KIT (FIVE PACK) 38015**

**ERGONOMIC
LIGHTWEIGHT
HANDLE** provides
better grip and
control.

-  **FINE SERRATION** on the blade is optimal to cut the periodontal ligament.
-  **THIN TAPERED** blade slides in to cut the PDL and reduces the breakage of the blade by reducing leverage on the overall blade.
-  **FLEXIBLE** blade takes the shape of the tooth to be able cut the PDL efficiently.



SURGICAL



38062 MOLT 2/4



37001 SINUS PLUGGER # 7

Visit our website
tbsinox.com
or call us
+1-844-TBS-INOX (827-4669)
to learn more!

FOLLOW US



/TBSINSTRUMENTS

38066 CASTRO X
NEEDLE HOLDER



Feature a twist hinge that allows your suture to just slide through without ever getting stuck in the hinge.



Castro X tips are crafted with **TUNGSTEN CARBIDE** to ensure optimal grip when grabbing a suture and needle for optimal grip.

TBS INOX USA

100% CRAFTED WITH
GERMAN STEEL

SURGICAL



MINNESOTA RETRACTOR 26002A

1-15 MARKING



IMPLANT PROBE 36001B

CONE SOCKET
WITH ENDO &
IMPLAN CALIPER



**TAWIL
MIRROR HANDLE 36010**

STRAIGHT,
TUGSTEN,
CARBIDE
TIPS



CASTRO X SCISSORS 38065



ASPIRATING FRINGE 38070

PERFORATED LARGE
END FOR BETTER
SUCTION



**PERIOSTEAL
ELEVATOR 38071**

SURGICAL



38072

**PRICHARD TISSUE
RETRACTOR**


WIDE AND NARROW ENDS
TO ELEVATE & RETRACT
TISSUE



38073

**PACKER PLUGGER
SERRATED**

2.5mm & 2mm
SERRATED ENDS
FOR PACKING BONE



38074

PERRY FORCEPS

TUNGSTEN CARBIDE
TIPS



38075

**ADSON TISSUE
FORCEPS 1x2**



38076

**MAYO NEEDLE
HOLDER**

TUNGSTEN CARBIDE
TIPS



38077

IRIS SCISSOR

TUNGSTEN
CARBIDE
TIPS

FOLLOW US



/TBSINSTRUMENTS

**TBS INOX
USA**

100% CRAFTED WITH
**GERMAN
STEEL**

SURGICAL

38500

TBS SURGICAL KIT W/CASSETTE

- 26002A 1 X MINNESOTA RETRACTOR
- 36001B 1 X IMPLANT PROBE
- 36010 1 X TAWIL MIRROR HANDLE
- 36020 1 X SURGICAL BOWL
- 37001 1 X SINUS PLUGGER
- 38002 1 X SMALL SERRATED CURETTE
- 38004 1 X MICRO SERRATED CURETTE
- 38018 1 X ANGLED PERIOTOME
- 38021 1 X STRAIGHT ELVATOME
- 38062 1 X MOLT 2/4
- 38063 1 X PUSH SCALPEL HANDLE
- 38065 1 X CASTRO X SCISSORS
- 38066 1 X CASTRO X NEEDLE HOLDER
- 38070 1 X ASPIRATING SYRINGE
- 38071 1 X PERIOSTEAL ELEVATOR
- 38072 1 X PRICHARD TISSUE RETRACTOR
- 38073 1 X PACKER PLUGGER
- 38074 1 X PERRY FORCEPS
- 38075 1 X ADSON TISSUE FORCEPS
- 38066 1 X MAYO NEEDLE HOLDER
- 38077 1 X IRIS SCISSORS

COMPLETE SURGICAL SETUP OF 21 INSTRUMENTS IN OUR EXCLUSIVE GOLD STERILIZABLE CASSETTE. TBS SURGICAL KIT INCLUDES 21 UNIQUELY DESIGNED INSTRUMENTS FOR ANY SURGICAL PROCEDURE. ALL 21 INSTRUMENTS FIT IN OUR EXCLUSIVE GOLD STERILIZABLE CASSETTE. TBS SURGICAL KIT IS THE MOST CONVENIENT AND EFFECTIVE WAY TO STORE, MANAGE AND STERILIZE YOUR SURGICAL INSTRUMENTS!

INSTRUMENT CARE & MAINTENANCE

1

FOLLOWING CLINICAL USE, ALL INSTRUMENTS MUST BE THOROUGHLY RINSED WITH WATER TO REMOVE ALL RESIDUE PRIOR TO INSTRUMENT PROCESSING.

NOTE: ABRASIVE BRUSHES OR MATERIALS SHOULD NOT BE USED TO CLEAN INSTRUMENTS.

2

INSTRUMENTS SHOULD BE PLACED INTO AN ULTRASONIC BATH FOR A PERIOD OF 20-30 MINUTES ACCORDING TO MANUFACTURER GUIDELINES TO ENSURE COMPLETE INSTRUMENT CLEANING.

NOTE: ONLY DEMINERALIZED WATER SHOULD BE USED IN AUTOCLAVE EQUIPMENT ACCORDING TO AUTOCLAVE MANUFACTURER INSTRUCTIONS.

3

ALL INSTRUMENTS MUST BE RINSED WITH WATER AND COMPLETELY DRIED BEFORE AUTOCLAVING. ALL STAINLESS-STEEL INSTRUMENTS SHOULD BE PLACED IN SEPARATE STERILIZATION POUCHES AND NOT TO BE MIXED WITH ANY CARBIDE INSTRUMENTS DURING INSTRUMENT PROCESSING.

FRINGS® CARE & MAINTENANCE


ALL FRINGS® MUST BE LUBRICATED WITH SURGICAL MILK TWICE A MONTH TO ENSURE OPTIMAL PERFORMANCE. WE RECOMMEND A SURGICAL MILK BATH EVERY TWO WEEKS AS PART OF YOUR INSTRUMENT MANAGEMENT PROCESS.

TITANIUM NITRIDE COATED INSTRUMENTS

TITANIUM COATINGS ARE USED TO INCREASE THE SURFACE HARDNESS OF INSTRUMENT TIPS AND REDUCE "PULLBACK" WHEN MANIPULATING COMPOSITE MATERIALS. PROPER CARE SHOULD BE TAKEN TO REMOVE ANY RESIDUAL COMPOSITE MATERIALS FROM THE BLADE WITHIN 5 MINUTES AFTER USE BY WIPING THE BLADE WITH ALCOHOL ON A 2x2 GAUZE. TITANIUM NITRIDE COATED INSTRUMENTS CAN BE CLEANED AND STERILIZED USING THE SAME RECOGNIZED ACCEPTABLE METHODS AS STAINLESS STEEL INSTRUMENTS

NOTE: NEVER USE ABRASIVES TO CLEAN THE SURFACE OF ANY TITANIUM COATED INSTRUMENTS.

**TBS INOX®
USA**



#findyourFRINGS

